



---







STAPLING  
RANGE

---

# CIRCULO 3R<sup>TM</sup>

- € Higher suture performance and anastomosis - 3 rows of staples
- € Single head design - outer ring 0.2 mm of the remaining staples rows - avoids excessive tissue tightening
- € Patented staple formation - offer a stronger suture, increasing tension, elimination mostly of the manual suturing reinforcements
- € Tissue ideal compression indicator - accurate and safe staple formation for variable tissue thickness
- € Black rubber insertions - extragrip and comfort during surgery
- € Button grip/release with wings - more control and accuracy during opening and closing of the anvil
- € Security system and security tab - ensures accurate formation in the various tissue thicknesses
- € Clear and intuitive marking

# REFERENCES

- CIR 21     CIRCULO 3R size 21
- CIR 24     CIRCULO 3R size 24
- CIR 25     CIRCULO 3R size 25
- CIR 26     CIRCULO 3R size 26
- CIRTD 24     CIRCULO 3R size 24 Tilt Dished
- CIRTD 26     CIRCULO 3R size 26 Tilt Dished
- 



# PRO-H VIEW™

€ Reinforced body - safer procedure  
Dark grey color - avoid reflections

€ Transparent cartridge - allow a better visual feedback

€ 0,75 mm suture clearance - maximizes hemostasis

€ Stainless steel anvil and stainless steel guide shaft - increase stability on the closing and firing procedures

€ High Volume housing - allow more tissue in the cartridge

€ Black rubber insertions - extragrip and comfort during surgery

€ Button grip/release with wings - more control and accuracy during opening and closing of the anvil

€ Security system and security tab - ensures accurate formation in the various tissue thicknesses





# REFERENCES

HVIEW-34  PRO-H-View size 34



# PRO-H 3R<sup>TM</sup>

- Ⓔ Higher suture performance and anastomosis - 3 rows of staples
- Ⓔ Single head design - outer ring 0.2 mm of the remaining staples rows - avoids excessive tissue tightening
- Ⓔ 3R Technology dramatically reduces bleeding during surgery but also post-operative bleeding
- Ⓔ Tissue ideal compression indicator - accurate and safe staple formation for variable tissue thickness
- Ⓔ Black rubber insertions - extragrip and comfort during surgery
- Ⓔ Button grip/release with wings - more control and accuracy during opening and closing of the anvil
- Ⓔ Security system and security tab - ensures accurate formation in the various tissue thicknesses
- Ⓔ 4 suture channels giving the possibility of a more accurate placement of the pull sutures
- Ⓔ Reinforced body - safer procedures  
Dark grey color - avoid reflections

# REFERENCES

H3R-34  PRO-H 3R size 34

H3R-36  PRO-H 3R size 36



# PRO-X<sup>TM</sup>

€ Retaining pin actuator

€ Reinforced body - safer procedure  
Dark grey color - avoid reflections

€ Orange release/safety button

€ Double trigger system - more intuitive, lock and fire system

€ Intuitive clamping/suture system of two triggers - (closure trigger/firing trigger) surface texturized to better grip

€ High quality stainless steel anvil, large jaw opening - accurate staple formation and regular compression

# REFERENCES

X03038S  PRO-X 30-3.8mm

X03045S  PRO-X 30-4.5mm

X06038S  PRO-X 60-3.8mm

X06045S  PRO-X 60-4.5mm

X09038S  PRO-X 90-3.8mm

X09045S  PRO-X 90-4.5mm

---

X03038R  RELOAD PRO-X 30-3.8mm

X03045R  RELOAD PRO-X 30-4.5mm

X06038R  RELOAD PRO-X 60-3.8mm

X06045R  RELOAD PRO-X 60-4.5mm

X09038R  RELOAD PRO-X 90-3.8mm

X09045R  RELOAD PRO-X 90-4.5mm



# PRO-X 3R<sup>TM</sup>

- € Higher suture performance and anastomosis - 3 rows of staples
- € Patented staple formation - higher suture tension
- € Retaining pin actuator
- € Intuitive clamping/suture system of two triggers - (closure trigger/firing trigger) surface texturized to better grip
- € Orange release/safety button
- € Reinforced body - safer procedures  
Dark grey color - avoid reflections
- € Ergonomic handle with a texturized skin to improve grip



# REFERENCES

X03025S  PRO-X 3R 30-2.5mm

X04545S  PRO-X 3R 45-45mm

X06045S  PRO-X 3R 60-4.5mm

X09045S  PRO-X 3R 90-4.5mm

---

X03025R  RELOAD PRO-X 3R 30-2.5mm

X04545R  RELOAD PRO-X 3R 45-45mm

X06045R  RELOAD PRO-X 3R 60-4.5mm

X09045R  RELOAD PRO-X 3R 90-4.5mm

---





# CORTE™

€ Pinless anvil for atraumatic placement

€ Safety lock-out - prevents accidental misfires before clamping procedure

€ Indicator of staple size in bright orange - prevents accidental placement of different staple sizes

€ Quick release button - easier stapler opening and reload operation

€ Black rubber inlays - extra grip and comfort

€ Clear graduations and indication of the cutting zone on the anvil

€ Reinforced body - safer procedure  
Dark grey color - avoid reflections

€ Adjustable orange firing button - can be operated by both sides

# REFERENCES

COR06038S  CORTE 60-38mm

COR06045S  CORTE 604.5mm

COR08038S  CORTE 80-3.8mm

COR08045S  CORTE 80-4.5mm

COR10038S  CORTE 100-3.8mm

COR10045S  CORTE 100-4.5mm

---

COR06038R  RELOAD CORTE 60-38mm

COR06045R  RELOAD CORTE 604.5mm

COR08038R  RELOAD CORTE 80-3.8mm

COR08045R  RELOAD CORTE 80-4.5mm

COR10038R  RELOAD CORTE 100-3.8mm

COR10045R  RELOAD CORTE 100-4.5mm



# CORTE 2<sup>TM</sup>

- € Pinless anvil for atraumatic placement
- € Safety lock-out - prevents accidental misfires before clamping procedure
- € Reinforced body - safer procedure  
Dark grey color - avoid reflections
- € Indicator of staple size in bright orange - prevents accidental placement of different staple sizes
- € Orange firing button
- € Clear and intuitive marking
- € Release handle - easier stapler opening and reload operation
- € Texturized grips - extra grip and comfort during surgery

# REFERENCES

COR206038S  CORTE 2 60-38mm

COR206045S  CORTE 2 60-4.5mm

COR208038S  CORTE 2 80-3.8mm

COR208045S  CORTE 2 80-4.5mm

COR210038S  CORTE 2 100-3.8mm

COR210045S  CORTE 2 100-4.5mm

---

COR206038R  RELOAD CORTE 2 60-38mm


COR206045R  RELOAD CORTE 2 60-4.5mm

COR208038R  RELOAD CORTE 2 80-3.8mm

COR208045R  RELOAD CORTE 2 80-4.5mm

COR210038R  RELOAD CORTE 2 100-3.8mm

COR210045R  RELOAD CORTE 2 100-4.5mm



# CURVO<sup>TM</sup>

€ Unique curved shape - ideal for procedures in the lower pelvic area

€ Retaining pin actuator - orange color

€ Ergonomic handle with a texturized skin for better grip

€ Orange release/safety button

€ Intuitive clamping/suture system of two triggers, (closure trigger/firing trigger) - surface texturized to better grip

€ Reinforced body - safer procedure  
Dark grey color - avoid reflections

# REFERENCES

CV03048S  CURVO 30-4.8mm

CV04048S  CURVO 40-4.8mm

---

CV03048R  RELOAD CURVO 30-4.8mm

CV04048R  RELOAD CURVO 40-4.8mm

---



# PELE<sup>TM</sup>

- € Simple and reasonable design - easy to operate
- € High sewing up and stapling speed during operation
- € Much less scar after cure
- € High quality stainless steel - avoids cross infection
- € Painless and convenient when removing the staples



# REFERENCES


PELEEXT  STAPLER REMOVER

PELE35W  SKIN STAPLER

---

---





Rua de Santo António  
30 A-Dos-Ralhados, 2710-006 Sintra  
LISBOA-PORTUGAL

